ODOFIN™


Sniffin' Sticks

Parosmie-Test - SSParot

burghart
messtechnik

DEUTSCH	3
ENGLISH	17
Legende/Legend	32

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Uetersener Straße 6
 D - 25488 Holm
 ☎ +49 4103 80076-0
 ☎ + 49 4103 80076-29
 ✉ info@burghart-mt.de
 🌐 www.burghart-mt.de

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ENGLISH**Content**

Intended use	18
User	18
Age group	18
Warnings	18
Contraindication.....	19
Patient information on the smell test	19
Information for users	20
Information for the Examiner	21
Requirement.....	21
Test procedure.....	21
Evaluation.....	22
Shelf Life.....	22
Storage	22
Disposal.....	23
Ingredients.....	24
Literature	25
Normative data	26

IMPORTANT
READ CAREFULLY BEFORE USE
KEEP FOR FUTURE USE

Intended use

The Sniffin' Sticks Parosmia test (SSParoT) measures the hedonic perception based on pairwise presented pleasant and unpleasant odors. The SSParoT indicates the presence of a qualitative olfactory impairment (Parosmia). The test can also monitor intra-individual changes of the hedonic perception, e.g. distorted olfaction during the course of a migraine. Moreover, the SSParoT may indicate a reduced hedonic perception that may occur due to neurological diseases (e.g. Morbus Parkinson).

User

The test is to be performed by healthcare professionals only.

Age group

The test is validated for people between age 18 to 35. Further normative data are needed for different age groups to allow an individual interpretation of results. A residual ability to smell is necessary to conduct the test.

Warnings

The olfactory sticks must not touch the patient's skin during the olfactory presentation. However, if the stick is inadvertently touched, it must be disposed of for hygienic reasons and because of the reduced shelf life of the odours and replaced by a new one.

All substances used in this test are not health hazardous in their used concentrations. Rarely allergic reactions show up after inhaling and skin contact. In case of allergic reaction contact the attending doctor and the manufacturer. Reevaluate!

The user shall wear odour-free gloves during the procedure.

Use of the sticks only in an odourless or well-ventilated environment.

Only use the stick for recommended procedures and not for writing purposes.

Caution: Ask for intolerances before each test.

Store the product safely and out of the reach of patients and children.

If the tip of the Sniffin' Sticks accidentally comes in contact with the gloves, the gloves need to be replaced to avoid contamination of other odours.

Contraindication


An allergy against one of the odorants used in the tests.

Patient information on the smell test


The patient shall be informed in the following manner: "Your sense of smell shall be tested with "Sniffin' Sticks" which contain different odors. Please smell each of the subsequently presented sticks and rate the odor on a scale where -4 is very unpleasant, 0 is neither pleasant nor unpleasant, and +4 is very pleasant. You also will be asked how intensive you perceive the odor (0, no perception, to 10, very intensive perception). The full test takes about 12 minutes."

Information for users

With the Parosmia test the patient rates the hedonic quality of 11 odor pairs on a 9-point scale from -4 to +4. Each odor pair is hedonically opposite, meaning the red stick contains a pleasant odor (E1) and the green stick an unpleasant odor (E2).

The hedonic range ($HR = E1 - E2$) and the hedonic direction ($HD = (E1 + E2) / 2$) of the pairs are calculated for evaluation. The hedonic range (HR, between -8 and +8) represents how different the opposite odors are perceived. The hedonic direction (HD, between -4 and +4) indicates if daily life odors are generally perceived pleasant or unpleasant. E.g., if both odors of a pair are perceived as very unpleasant, the hedonic range (HR) is small and the hedonic direction (HD) is negative, which might indicate a parosmia.

The intensity should be larger than 3 on the numerical rating scale to attain significant results. Otherwise, the hedonic quality is mainly rated 0.

The Parosmia test may only be used in connection with the user manual, the holder and the Odofin "Sniffin' Sticks". For refilling, only refill sets from the manufacturer are allowed to be used.

Information for the Examiner

Requirement

The following criteria must be taken into account:

- To obtain the most reliable results possible, the patient should consume nothing other than water at least 30 minutes before the start of the test.
- Chewing gum, sweets or cigarettes must also be avoided.
- During the test, the patient is not to be given any indication of the quality of the odor. There is no correct or wrong answer because it is asked for the individual and subjective perception.
- The Identification test may only be used in combination with instructions for use, holders and sticks. Only the refill set supplied by the manufacturer may be used as a refill set.

Test procedure

Before presenting the Sniffin' Stick, hand out the corresponding multiple-choice card. Pull the Sniffin' Stick out of the holder and remove the cap carefully. Sometimes a gentle rotation of the cap is necessary.

For testing, the opened odour stick is held 2 cm centred in front of both nostrils and the patient is asked to smell, e.g., using the sentence "Please sniff now". Each stick should be presented for 3-4 seconds. The interval between the different presentations should be around 30 seconds.

Ask the patient to rate the odor intensity on the numerical rating scale (NRS) and subsequently determine the Hedonic Quality on the 9-point scale. Write down the answers on the evaluation sheet (red stick E1 and green stick E2). Subsequently present the 11 pairs to the patient with an interval of around 30 seconds.

Please put the individual stick back into the block directly after use and with the cap closed. The sticks must be stored in an upright "overhead" position (white cap downwards). This prevents release of fragrance into the room.

Evaluation

For every odor pair the hedonic range (HR) and hedonic direction (HD) must be calculated from the answers on the evaluation sheet. If the odor intensity was rated too low (<3) the odor pair does not count for further analysis. The mean value of the remaining HRs and HDs must be calculated and compared to the normative values in table 1 (HR) and table 2 (HD).

Shelf Life

The sticks are labelled with a lot-number and a best before date. Use beyond the expiry date may lead to incorrect results.

Storage

Refer to the symbols for how to store the test (dry, protect from sunlight). The temperature should be between 10°C and 30°C.

The sticks must be stored in an upright "overhead" position (white cap downwards).

Caution:

Storage in the fridge can lead to a change in the intensity of the odour.

Disposal

Unusable sticks must be disposed of in their sealed state in the household waste.

Ingredients**Red Sniffin' Sticks (hedonic quality: pleasant)**

Pen	Fragrance
No. 1	Peppermint
No. 2	Apple
No. 3	Pineapple
No. 4	Banana
No. 5	Peach
No. 6	Coconut
No. 7	Caramel
No. 8	Raspberry
No. 9	Ice bonbon
No. 10	Lemon
No. 11	Orange
Solvent	Diethyl phthalate
Solvent	Propylenglycolum

Green Sniffin' Sticks (hedonic quality: unpleasant)

Pen	Fragrance
No. 1	Fish
No. 2	Garlic
No. 3	Turpentine
No. 4	Clove
No. 5	Butter
No. 6	n-butyric acid solution
No. 7	Iso-butyric acid solution
No. 8	Indol solution
No. 9	Skatole solution
No. 10	Civet solution
No. 11	Valeric acid solution
Solvent	Diethyl phthalate
Solvent	Propylenglycolum

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Assessment of odor hedonic perception: the Sniffin' sticks parosmia test (SSParoT)

Sci Rep. 2020 Oct 22;10(1):18019.

doi: 10.1038/s41598-020-74967-0.

PMID: 33093474; PMCID: PMC7581750.

Helpfenstein C, Strupf M, Stefke A, Fraunberger B, Renner B, Suchantke I, Rothermel M, Messlinger K, DeCol R, Namer B.

Cyclic changes of sensory parameters in migraine patients

Cephalalgia. 2022 Oct;42(11-12):1148-1159.

doi: 10.1177/03331024221097932.

PMID: 35514204.

Mrochen A, Marxreiter F, Kohl Z, Schlachetzki J, Renner B, Schenk T, Winkler J, Klucken J. **From sweet to sweat:**

Hedonic olfactory range is impaired in Parkinson's disease

Parkinsonism Relat Disord. 2016 Jan;22:9-14.

doi: 10.1016/j.parkreldis.2015.09.035.

PMID: 26627940.

Normative data

Following tables contain the normative data for HR and HD compiled from Liu et al. 2020 (see “Literature”). Study participants were between 18 and 35 years old. Further normative data are needed for different age groups to allow an individual interpretation of results.

Pair	Hedonic Range		10. Perzentile	
	♂	♀	♂	♀
1	5.6	5.5	2.0	3.0
2	4.9	4.6	1.4	2.0
3	3.8	3.5	0.0	0.0
4	3.8	2.1	1.0	0.0
5	4.0	5.0	1.0	4.0
6	4.2	4.0	2.0	0.0
7	3.7	4.2	1.0	1.2
8	3.1	3.9	-1.0	0.8
9	4.1	4.7	1.0	1.8
10	3.5	4.3	0.0	0.8
11	4.8	4.7	2.0	3.0
Ø 1-11	4.1	4.4	1.0	2.0

Table 1: Normative values for the hedonic range (HR)

The average hedonic range (HR) is 4.1 for males and 4.4 for females. The 10th percentile is 1.0 for males and 2.0 for females, meaning that only 10% of normosmic people perceived the hedonic range between pairs as very small.

Example A:

The average HR of a male patient is 1.5. He regarded most odors as intense enough (>3). The value is above the 10th percentile of 1.0, therefore there is no indication of a parosmia.

Example B:

The average HR of a female patient is 1.5. She regarded most odors as intense enough (>3). The value is below the 10th percentile of 2.0, therefore there is an indication for a parosmia.

Pair	Hedonic Direction		10. percentile	
	♂	♀	♂	♀
1	-1.3	-1.4	-2.0	-2.0
2	-0.6	-0.9	-1.5	-2.5
3	0.2	0.0	-1.0	-2.0
4	1.4	2.0	-0.5	0.5
5	0.5	-0.6	-0.5	-0.6
6	0.2	-0.6	-1.0	-0.6
7	0.0	-1.2	-1.5	-1.2
8	-0.2	-1.5	-1.5	-1.5
9	-0.1	-2.0	-1.5	-2.0
10	-0.3	-2.0	-1.5	-2.0
11	-0.6	-1.6	-1.5	-1.6
Ø 1-11	-0.1	-0.1	-1.5	-2.0

Table 2: Normative values for the hedonic direction (HD)

The average hedonic direction (HD) is -0.1 for males and females, so about 0. The 10th percentile for males is -1.5 and for females -2.0, meaning that only 10% of normosmic people regard the mean hedonic quality of the pairs as relatively unpleasant.



Example A:

The average HD of a male patient is -1.9. He regarded most odors as intense enough (>3). The value is below the 10th percentile of -1.5, therefore there is an indication for a parosmia.

Example B:

The average HD of a female patient is -1.9. She regarded most odors as intense enough (>3). The value is just over the 10th percentile of -2.0, therefore there is no indication for a parosmia.


**Copy template for the numerical rating scale (NRS)
and the 9-point hedonic-scale:**

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